

CERTIFICATE OF ANALYSIS

Product Name : Propylene Glycol

Lot code : 2415

MF : C₃H₈O₂

Mfg. Date : May 2026

Analyzed On : 28-April-2026


MW : 76.10

Exp. Date : June 2031

CAS No. : 57-55-6

Tests	Specifications	Results
Appearance	Clear colourless liquid	Clear colorless liquid
Purity (by GC)	Min 99%	99.80%
Density (d _{20°C/4°C})	1.035 - 1.040	1.037
Refractive index (20°C; 589 nm)	1.4320 - 1.4330	1.4323
10% aqueous solution	Clear	Clear
Water	Max 0.3%	0.14%
Acidity or alkalinity	Max 0.3 ml N%	0.05 ML N%

This above product complies as per the specifications of : **MOLATOM SCIENTIFIC. LTD.**



MOLATOM SCIENTIFIC LTD.

Website: [www.https://www.molatomscientific.com](https://www.molatomscientific.com)