

CERTIFICATE OF ANALYSIS

Product Name : TRIETHYLENE GLYCOL 98% For Synthesis

Lot code : 3120

Mol. Formula : C₆H₁₄O₄

Code No. : 00347

Mfg. Date : Jan-2026

Release Date : 08-Jan-2026


Mol. Weight : 150.18

CAS No. : 112-27-6

Exp. Date : Dec-2030

Tests	Specifications	Results
Appearance	Clear colorless to light yellow liquid	light yellow liquid
Purity (by GC)	Min 98%	99.0%
Density (d _{20°C/4°C})	1.123 - 1.124	1.123
Water	Max 0.3%	0.18%
Refractive index (20°C; 589 nm)	~1.4559	1.4542

This above product complies as per the specifications of : **MOLATOM SCIENTIFIC. LTD.**



MOLATOM SCIENTIFIC LTD.

Website: [www.https://www.molatomscientific.com](https://www.molatomscientific.com)